


VACCINE TRANSPORT COOLERS MUST BE PRE-CONDITIONED BETWEEN $+2.0^{\circ}\text{C}$ to $+8.0^{\circ}\text{C}$ PRIOR TO VACCINE PICKUP

Complete these 3 Easy Steps:

- 1) Arrange vaccine transport cooler packing materials and wait approximately 30 minutes for temperatures to stabilize between $+2.0^{\circ}\text{C}$ to $+8.0^{\circ}\text{C}$.
- 2) Clear the thermometer's memory.
- 3) Proceed to your vaccine pickup location. York Region staff will verify that the vaccine transport cooler is within $+2.0^{\circ}\text{C}$ to $+8.0^{\circ}\text{C}$ before releasing your vaccine order.



Reminder: Check your cooler's temperature upon arrival at your office. Immediately contact the Vaccine Inventory program at 1-877-464-9675 ext. 74065 if the temperature is below $+2.0^{\circ}\text{C}$ or above $+8.0^{\circ}\text{C}$.

22-5148