

York Region Infection Prevention And Control Lapse Report

Initial Report				-	
Premise/Facility under investigation (name and address Sil-Dom Dental Clinic Unit 7-10555 Jane Street)				
Maple, Ontario L6A 3L1					
Type of Premises/Facility Dental Clinic					
Date Board of Health became aware of IPAC lapse (yyyy/mm/dd) 2019/10/2			Date of Initial Report posting (yyyy/mm/dd) 2019/10/07		
Date of Initial Report update(s) (if applicable) (yyyy/mm/dd)			How the IPAC lapse was identified Complaint		
Summary Description of the IPAC Lapse During an inspection on October 2, 2019 York Region Public Health noted the following: dental equipment/devices were not disassembled and in the open position for sterilization; dental equipment/devices were bound together during ultrasonic cleaning; sterilizer Universal Serial Bus (USB) record was not checked to ensure sterilization had been achieved for each sterilized load; sterilizer was not tested with a biological indicator (BI) in a process challenge device (PCD) for each type of cycle and for each day the sterilizer was used; a Type 5 or 6 chemical indicator in a PCD was not used to justify the release of sterilized loads of dental equipment/devices prior to BI results.					
IPAC Lapse Investigation	Yes	No	N/A	Please provide further details/steps	
Did the IPAC lapse involve a member of a regulatory college?	X				
If yes, was the issue referred to the regulatory college?	X			Royal College of Dental Surgeons of Ontario, College of Dental Hygienists of Ontario	
Were any corrective measures recommended and/or implemented?	X				
Please provide further details/steps	At the time of the inspection, operator was educated on correct infection prevention and control (IPAC) practices. Operator is required to ensure that: Dental equipment/devices are disassembled and in the opened position during the cleaning and sterilization process unless otherwise stated in the manufacturer's instructions Dental equipment/devices during cleaning are separated and not bound together in the ultrasonic cleaning solution Sterilizer mechanical display, print out or Universal Serial Bus (USB) record is checked to ensure sterilization has been achieved and verify this by initialing in a log book for each sterilized load Sterilizer is tested with a biological indicator (BI) in a process challenge device (PCD) each day the sterilizer is used and with each type of cycle used that day Dental equipment/devices are only released when the BI results are available; if quarantine pending BI results is not possible, evaluation of a Type 5 or 6 chemical indicator (CI) in a PCD and the specific cycle parameters are used to justify the release of sterilized loads The operator demonstrated the correct IPAC practices for the above items after education was provided.				
Date any order(s) or directive(s) were issued to the owner/operator (if applicable) (yyyy/mm/dd) 2019/10/2					
Initial Report Comments and Contact Information					
Any additional Comments: (Please do not include any personal information or personal health information)					



York Region Infection Prevention And Control Lapse Report

If you have any further questions, please contact
Health Connection

Telephone Number

1-800-361-5653

Email Address
Health.inspectors@york.ca

Final Report

Date of Final Report posting (yyyy/mm/dd)

2019/10/07

Date any order(s) or directive(s) were issued to the owner/operator (if applicable) (yyyy/mm/dd) 2019/10/02

Brief description of corrective measures taken

An inspection report was provided to the clinic on Oct. 2, 2019 noting corrective measures to be implemented. An Order under the Health Protection and Promotion Act was issued outlining corrective measures to be implemented. The operator was educated on-site. A re-inspection was conducted on Oct. 4, 2019 where it was noted by York Region Public Health that all corrective measures have been implemented.

Date of all corrective measures were confirmed to have been completed (yyyy/mm/dd) 2019/10/4

Final Report Comments and Contact Information

Any Additional Comments: (Please do not include any personal information or personal health information)

If you have any further questions, please contact Health Connection

Telephone Number Email Address

1-800-361-5653 Health.inspectors@york.ca