

PROPER CLEANING AND DISINFECTION PRACTICES

Cleaning

- Cleaning must always be the first step to remove dirt and debris from a surface and is necessary for a disinfectant to be effective
- Clean with a detergent, water, and friction and clean from least contaminated to most contaminated areas

Disinfectants

- Disinfectants are applied to a clean surface in order to kill disease-causing germs
- Disinfectants must have a drug identification number (DIN) if approved for use in Canada (common household bleach and isopropyl alcohol are the only exceptions)
- Always follow manufacturer's instructions for use (MIFU). Read label for direction on: dilution and mixing, personal protective equipment (PPE) needed (e.g., gloves, goggles), surfaces appropriate for use, contact time, efficacy on specific organisms, and rinsing requirements
- There are a variety of disinfectants in the market. Choose a disinfectant that is compatible with your surfaces and with contact times that fit your needs
- Check the expiry date. If a product has expired, do not use. Discard expired product safely or return to manufacturer
- Ensure the concentration of disinfectant is correct before use (i.e. use test strips)
- Toys that will be mouthed should be rinsed thoroughly with water following disinfection
- Do not use antiseptic wipes and other products intended for skin (i.e. alcohol-based hand rubs) on surfaces

Cleaning and disinfecting wipes

- Follow manufacturer's recommendations
- Wipes may become dry (improper storage or during use) due to fast drying properties before contact time is achieved
- Wipes are not recommended as a routine cleaning/ disinfectant tool, especially for heavily soiled surfaces
- They can be used for items that cannot be soaked and for small items that must be disinfected between uses
- Ensure the surface or item remains wet with the product for the required contact time (additional wipes may be needed)
- Wipes must be kept wet and should be discarded if they become dry

Blood and body fluid spills

- Wipe spills immediately- use disposable towels to remove most of the organic matter, clean the area and then disinfect the spill area
- See the chart below for examples of disinfectants to use depending on volume of blood/ body fluid spill

When to Clean and Disinfect	Examples of Active Ingredients/ Disinfectant Products	Contact Time (minutes)	Where to Clean and Disinfect
1. Everyday use (non-outbreak) 2. Minor blood/ body fluid spill (drops of fluid) Effective against: <i>Vegetative bacteria and enveloped viruses</i> <ul style="list-style-type: none"> • Staphylococcus aureus (includes MRSA) • Streptococcus • Salmonella • Vancomycin Resistant Enterococcus (VRE) • Human Immunodeficiency Virus (HIV) • Respiratory Syncytial Virus (RSV) • Influenza Virus • Pseudomonas Aeruginosa • Herpes 	100 ppm bleach solution <i>(Everyday use; non-outbreak)</i>	10	Surfaces: <ul style="list-style-type: none"> • Door knobs • Hand rails • Chairs • Tables • Elevator buttons • Telephones • Counter tops • Sink faucet handles • Toys • Commode chairs • Shared play equipment • Vinyl mattress covers • Floor mats • Water fountains • Diaper change stations Equipment: <ul style="list-style-type: none"> • Blood pressure cuffs • Thermometers • Stethoscope
	1,000 ppm bleach solution <i>(Minor blood/ body fluid spill)</i>	10	
	Quaternary Ammonium Compounds (QUATS) <i>(i.e. Lysol®, ED- Everyday Disinfectant, Quato 78 Plus™, A-3®, Swish Clean and Green™)</i>	MIFU	
	3% Hydrogen Peroxide	10	
	70-90% Alcohol <i>(Ethyl or Isopropyl)</i> <i>(For soaking)</i>	10	
	Zochlor	MIFU	
	Benefect®	10	
1. Outbreak situation 2. Major blood/ body fluid spill 3. Confirmed viral or bacterial infection of pathogens listed below (non-outbreak situation) Effective against: <i>Mycobacteria, enveloped and non-enveloped viruses and fungi</i> <ul style="list-style-type: none"> • Mycobacteria tuberculosis • Norovirus • Hepatitis A Virus • Rotavirus • Coxsackie Virus/ Hand, Foot and Mouth Disease • Rhinovirus/ Common Cold • Candida 	5,000 ppm bleach solution <i>Also a sporicidal</i> <i>(see below for recipe)</i>	10	
	6% Hydrogen Peroxide	30	
	Enhanced Action Formulation Hydrogen Peroxide	MIFU	
	Zochlor	MIFU	

Note: York Region Community and Health Services does not endorse any of the examples of brand name products listed above.

Bleach (Sodium Hypochlorite) Solutions

- Use undiluted household bleach (5.25% or ~50,000 ppm) when making the solutions in the chart below
- When making bleach solutions, add bleach to water- do not add water to bleach
- Store bleach solutions in closed containers, away from heat and light
- Bleach solutions should be properly labelled
- Online dilution calculator available from Public Health Ontario at the following link:
<https://www.publichealthontario.ca/en/health-topics/environmental-occupational-health/water-quality/chlorine-dilution-calculator>

Parts per million (ppm) Concentration	Recipes (~ = approximately, t = teaspoon, tb = tablespoon)
100 ppm (1:500, 0.01%)	Mix 2 ml (0.4 t) of bleach with ~1 L (4 cups) of water
1,000 ppm (1:50, 0.1%)	Mix 20 ml (4 t) of bleach with ~1 L (4 cups) of water
5,000 ppm (1:10, 0.5%)	Mix 100 ml (6 ¾ tb) of bleach with ~1 L (4 cups) of water

Community and Health Services
Public Health

1-877-464-9675
TTY 1-866-512-6228
york.ca

